




F12x8.2 F12x8.2

12 x 8 Strips of 200 μ L tubes



Fixed-Angle Rotor, Aluminum

Major applications: Minimizing post-PCR* condensation, post-reaction clean-up.

Max. rpm	Max. g	Number of Tubes Volume/Size	Approximate Accel/Decel Time (min:sec)
14,000	20,800 (outer row) 19,590 (inner row)	12 x 8 Strips 200 μ L	0:38/0:35

For use in Microfuge® 22R microcentrifuge.

2

No. 369534. F12x8.2 Rotor Assembly.

Tubes

Tube Style/Material	Nominal Volume per Tube (μ L)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed
PCR Tubes	200	†	—	—	—	—	—	—	14,000

* PCR is covered by patents owned by F. Hoffman-La Roche, Inc.

† Commercially available.

Rotor Replacement Parts

369535 Rotor	Lid Assembly
368896 Rotor	O-ring
369537 Lid	O-ring
365806 T	ie-down Screw
365636 T	-handle Rotor Wrench