

LED LOUPES LIGHT SOURCE

The Miltex® portable dental and surgical LED light source produces amazing brightness in a completely portable light. The small power pack is lightweight and uses Lithium Ion cells that can be quickly and easily recharged. The portable LED can be attached to most loupes systems.

Features	Benefits
Highest intensity of any portable dental and surgical headlight (80,000 Lux @ 14 inches)	Provides pure white light for true color rendition
Power pack design with external battery cells for quick changing Li Ion batteries (1 included)	Provides 4.5 hours of use on a single charge
Variable light intensity control	Allows for professional to adjust light as needed
Portable Headlamp Weighs only 0.7 ounces	Lightweight reduces pressure
Quick charger - 90% charge in 90 minutes	Reduced charging time
Flip up dental composite filter available	Prevents composite materials from curing



- Each unit comes complete with 1 each of the following:
- 1. LED Headlight
 - 2. Battery
 - 3. Battery Charger
 - 4. Battery Power Pack (attaches to belt)
 - 5. AC Adapter
 - 6. Attach mount (varies per REF)