

Combo V24T Centrifuge

Accomplish both test tube and microhematocrit separations with the Combo V24T centrifuge. The Combo V24T spins test tubes, capillary tubes or microtubes by simply changing rotors — no tools required! Achieve microhematocrit separations in three minutes at 12,000 rpm; then accurately determine results with the included reader disc. Other features include digital speed and time display, and all-metal body construction.

ADVANTAGES

- One unit for all centrifuge needs
- Fast spin, fast stop, fast read
- Variable speed
- Digital controls
- Quiet vibration-free operation
- Spin Microhematocrit or Test Tubes
- Easy-Change Rotors

FEATURES

- Rubber feet to prevent slipping
- Impact-resistant safe-view lid
- All-metal body construction
- Cooling fan for motor compartment
- Heavy duty high-speed motor
- Safety shut off switch in lid
- Locking lid
- Wide mouth for easy clean-up
- Digital readout

SPECIFICATIONS

12,000rpm on High speed setting

- Microhematocrit rotor: 14,500g
- Microtube rotor: 12,900g

4,000rpm on Low speed setting

- Test tube rotor: 2,415g

110v power (220v available)

CONTROLS

- Auto brake for quick stop
- High/low speed settings
- Digital programmable speed control
- Digital programmable time control

WEIGHT AND DIMS

Height: 11" (277mm)
Length: 17" (432mm)
Width: 14" (355mm)
Weight: 31.3 lbs (14.2kg)



Combo V24T Centrifuge Complete Kits

Model #	Description
CMC-1224-77T1	V24T centrifuge with 12-place test tube and 24-place microhematocrit rotors, 110v
CMC-24AV-15T1	V24T centrifuge with 24-place all metal test tube rotor, 110v
CMC-1224-77T2	V24T centrifuge with 12-place test tube and 24-place microhematocrit rotors, 220v
CMC-24AV-15T2	V24T centrifuge with 24-place all metal test tube rotor, 220v

CUSTOM BUILD YOUR COMBO V24T CENTRIFUGE



CMP-BX7N-77T1 (110v), CMP-BX7N-77T2 (220v)
Combo V24T Centrifuge Unit (rotors not included)

Choose from the rotor(s) below:



CMP-RT24-75HE

Spin standard 75mm capillary tubes with the 24-place microhematocrit rotor. Results can be read directly from the rotor with the included reader disk.



CMP-RT24-02MT

The 24-place microtube rotor holds 1.5 to 2 ml sized tubes.



CMP-RT12-15TT

The 12-place test-tube rotor holds 3ml-15ml sized tubes with included tube shields and inserts for smaller tubes.



CMP-RT24-15TT

The 24-place test-tube rotor holds 3ml-15ml sized tubes with included tube shields and inserts for smaller tubes.