Water Bath Blankets

LW Scientific's Water Bath Blankets are an incredible addition to your lab or clinic and save you money and time. Use them instead of a lid on your water bath. Immerse your tubes through the ball thermal layer, and never have to worry about removing a lid and replacing it! The thermal ball layer retains the water vapor, and the heat! Because it is so efficient at retention, it lowers the work load on your water bath's heaters and motors, making them last longer and use less energy. LW Scientific's Bath Blankets are available in various sizes and quantities.



Water Bath Blankets

	Size	Weight	Wall Thickness	Diameter Tolerance	Surface Coverage
10mm	0.3937 in.	0.2 grams	1.2 mm	+/- 0.25 mm	1,076 pcs/ft ²
20mm	0.787 in.	0.95 grams	1.2 mm	+/- 0.25 mm	250 pcs/ft ²
38mm	1.5 in.	4 grams	1.2 mm	+/- 0.50 mm	70 pcs/ft ²
50mm	2 in.	8 grams	1.2 mm	+/- 0.25 mm	43 pcs/ft ²





P 770.270.1394 F 770.270.2389

865 Marathon Parkway Lawrenceville GA 30046

SLS-7.5.3-L-074 | Rev 1 LWScientific.com