

IDC/GAF



DESCRIPTION

Infectious disease control goggle. A complete foam cushion goggle assembly with elastic strap attachment. Designed for splash, smoke & particulate protection for the eyes.

SPECIFICATION

MATERIAL	Mylar (Anti-Fog both sides); Elastic; High performance foam
WEIGHT	9 g (0.33 oz)
COLOR	Clear
STANDARD	OSHA 1910.1030, EN 166

ORDER DATA

PRODUCT NAME	IDC/GAF
PRODUCT NUMBER	4630006

APPLICATION

Evacuation, EMT's, Hospital Staff, Labs, Dental, etc.

Declaration of Conformity:

[DoC Part 1: Cover Sheet](#)

[DoC Part 2: Attachment](#)